

Classic Style Stethoscope Air Tague Holder Assembly

Thank you for purchasing the Stethoscope Air tag holder See below for assembly:

Place Stethoscope tube in center channel if it fits snuggle the center insert in the first picture may not be needed.


Place air tag in the circular area in the base. The pins on the top align with the holes in the base see picture 2 yellow circles.

Place screw in the screw holes and tighten down.

If you have any issues let us know shamrockwood315@gmail.com.

www.Shamrockwood315.com

Picture 1



Picture 2

